
Centrifuge Tubes

CAT No.	TUBE DESCRIPTION
MQ-BS-15-25	15 ml Centrifuge Tube with flat Cap, Sterile, Bag of 25
MQ-BS-50-25	50 ml Centrifuge Tube with flat Cap, Sterile, Bag of 25





Description

The MOLEQULE-ON 15ml and 50ml centrifuge tubes are manufactured from FDA approved resin and free from RNase, DNase, DNA, Endotoxin and Pyrogen. The tubes are also free from heavy metals and additives. These tubes can be centrifuged upto 20,000 rcf and can be stored from -80°C to +121°C. The lids are specially designed to to prevent leakage. Frosted area is available for labeling. No sharp points or edges, ultra smooth design that does not catch gloves. Clear and precise graduation with immaculate visibility to perform all types of experiments. The tubes can be can be fit in all standard and high capacity centrifuge rotors. These are sterile thus suitable for any routine and sensitive laboratory and diagnostic work.